

Package ‘fmriqa’

May 8, 2026

Type Package

Title Functional MRI Quality Assurance Routines

Version 0.5.0

Date 2026-02-20

Description Methods for performing fMRI quality assurance (QA) measurements of test objects. Heavily based on the fBIRN procedures detailed by Friedman and Glover (2006) <[doi:10.1002/jmri.20583](https://doi.org/10.1002/jmri.20583)>.

BugReports <https://github.com/martin3141/fmriqa/issues>

License GPL-3

RoxygenNote 7.3.3

Imports viridisLite, RNifti, ggplot2, reshape2, gridExtra, grid, tidy, optparse, tcltk, RcppEigen, imager, pracma

Encoding UTF-8

Suggests testthat, covr, knitr, rmarkdown

VignetteBuilder knitr

NeedsCompilation no

Author Martin Wilson [cre, aut]

Maintainer Martin Wilson <martin.wilson111@gmail.com>

Repository CRAN

Date/Publication 2026-02-20 15:40:02 UTC

Contents

| | |
|----------------------------|---|
| fmriqa-package | 2 |
| combine_res_glob | 2 |
| run_fmriqa | 3 |
| run_fmriqa_glob | 4 |

| | |
|--------------|----------|
| Index | 6 |
|--------------|----------|

fmriqa-package

fmriqa: fMRI quality assurance routines

Description

The fmriqa package provides an implementation of the fMRI quality assurance analysis protocol detailed by Friedman and Glover (2006) <doi:10.1002/jmri.20583>.

Details

#' To learn more about fmriqa, start with the vignettes: `'browseVignettes(package = "fmriqa")'`

For a full list of functions: `'help(package=fmriqa, help_type="html")'`

Author(s)

Maintainer: Martin Wilson <martin.wilson111@gmail.com>

See Also

Useful links:

- Report bugs at <https://github.com/martin3141/fmriqa/issues>

combine_res_glob

Combine fmriqa csv result files into a dataframe

Description

Combine fmriqa csv result files into a dataframe

Usage

```
combine_res_glob(pattern)
```

Arguments

pattern glob pattern to match the csv result files (eg "*.csv")

Value

a dataframe of results

`run_fmriqa`*Run fMRI quality assurance procedure on a NIfTI data file*

Description

Run fMRI quality assurance procedure on a NIfTI data file

Usage

```
run_fmriqa(  
  data_file = NULL,  
  roi_width = 21,  
  slice_num = NULL,  
  skip = 2,  
  tr = NULL,  
  pix_dim = NULL,  
  poly_det_ord = 3,  
  spike_detect = FALSE,  
  x_pos = NULL,  
  y_pos = NULL,  
  plot_title = NULL,  
  last_vol = NULL,  
  gen_png = TRUE,  
  gen_res_csv = TRUE,  
  gen_pdf = FALSE,  
  gen_spec_csv = FALSE,  
  png_fname = NULL,  
  res_fname = NULL,  
  pdf_fname = NULL,  
  spec_fname = NULL,  
  verbose = TRUE,  
  bg_smooth = 12,  
  bg_shrink = 25,  
  t1_canny = NULL,  
  t2_canny = NULL  
)
```

Arguments

| | |
|------------------------|--|
| <code>data_file</code> | input data in nifti format, a file chooser will open if not set |
| <code>roi_width</code> | roi analysis region in pixels (default=21) |
| <code>slice_num</code> | slice number for analysis (default=middle slice) |
| <code>skip</code> | number of initial volumes to exclude from the analysis (default=2) |
| <code>tr</code> | override the TR detected from data (seconds) |
| <code>pix_dim</code> | override the x,y,z pixel dimensions (mm) detected from data eg pixdim=c(3,3,3) |

| | |
|--------------|--|
| poly_det_ord | polynomial order used for detrending (default=3) |
| spike_detect | generate k-space spike-detection plot (default=FALSE) |
| x_pos | x position of ROI (default=center of gravity) |
| y_pos | y position of ROI (default=center of gravity) |
| plot_title | add a title to the png and pdf plots |
| last_vol | last volume number to use in the analysis |
| gen_png | output png plot (default=TRUE) |
| gen_res_csv | output csv results (default=TRUE) |
| gen_pdf | output pdf plot (default=FALSE) |
| gen_spec_csv | output csv of spectral points (default=FALSE) |
| png_fname | png plot filename |
| res_fname | csv results filename |
| pdf_fname | pdf plot filename |
| spec_fname | csv spectral data filename |
| verbose | provide text output while running (default=TRUE) |
| bg_smooth | amount to smooth background image before calculating the maximum BG percent metric (default=12mm) |
| bg_shrink | amount to shrink the BG image away from the object to avoid residual object signal in the maximum BG percent metric (default=25mm) |
| t1_canny | threshold for weak edges for Canny edge detector (defaults to auto detect) |
| t2_canny | threshold for strong edges for Canny edge detector (defaults to auto detect) |

Value

dataframe of QA metrics

Examples

```
fname <- system.file("extdata", "qa_data.nii.gz", package = "fmriqa")
res <- run_fmriqa(data_file = fname, gen_png = FALSE, gen_res_csv = FALSE, tr = 3)
```

run_fmriqa_glob

Run fMRI quality assurance procedure on a set of NIfTI data files

Description

Run fMRI quality assurance procedure on a set of NIfTI data files

Usage

```
run_fmriqa_glob(pattern, ...)
```

run_fmriqa_glob

5

Arguments

| | |
|---------|---|
| pattern | glob expression to match analysis files |
| ... | options to pass to run_fmriqa function |

Index

`combine_res_glob`, [2](#)

`fmriqa (fmriqa-package)`, [2](#)

`fmriqa-package`, [2](#)

`run_fmriqa`, [3](#)

`run_fmriqa_glob`, [4](#)