

# Package ‘isoorbi’

May 8, 2026

**Type** Package

**Title** Process Orbitrap Isotopocule Data

**Version** 1.5.2

**URL** <https://isoorbi.isoverse.org/>, <https://github.com/isoverse/isoorbi>

**BugReports** <https://github.com/isoverse/isoorbi/issues>

**Depends** R (>= 4.4.0)

**Imports** utils (>= 4.1.0), stats (>= 4.1.0), methods, rlang (>= 1.0.0), cli (>= 3.6.0), glue, withr (>= 3.0.0), lifecycle (>= 1.0.0), tidyr (>= 1.2.0), tibble, tidyselect (>= 1.2.0), dplyr (>= 1.1.1), ggplot2 (>= 3.4.0), scales (>= 1.2.1), readr (>= 2.1.0), readxl, openxlsx, purrr, prettyunits (>= 1.2.0), arrow (>= 21.0.0.0), knitr (>= 1.5.0)

**Suggests** devtools, fansi, rmarkdown, testthat (>= 3.0.0), vdiff, forcats

**Description** Read and process isotopocule data from an Orbitrap Isotope Solutions mass spectrometer. Citation: Kantnerova et al. (Nature Protocols, 2024).

**License** MIT + file LICENSE

**Encoding** UTF-8

**VignetteBuilder** knitr

**RoxygenNote** 7.3.2

**Config/testthat/edition** 3

**NeedsCompilation** no

**Author** Caj Neubauer [aut, cre, cph] (ORCID:

<https://orcid.org/0000-0002-5348-5609>),

Sebastian Kopf [aut] (ORCID: <https://orcid.org/0000-0002-2044-0201>),

Kristýna Kantnerová [aut] (ORCID:

<https://orcid.org/0000-0001-6259-3225>)

**Maintainer** Caj Neubauer <caj.neubauer@colorado.edu>

**Repository** CRAN

**Date/Publication** 2025-10-03 21:10:02 UTC

## Contents

isoorbi-package . . . . .	3
orbi_add_blocks_to_plot . . . . .	3
orbi_adjust_block . . . . .	4
orbi_aggregate_raw . . . . .	6
orbi_analyze_shot_noise . . . . .	7
orbi_calculate_ions . . . . .	7
orbi_calculate_ratios . . . . .	8
orbi_calculate_summarized_ratio . . . . .	9
orbi_check_isoraw . . . . .	10
orbi_default_theme . . . . .	11
orbi_define_basepeak . . . . .	11
orbi_define_blocks_for_dual_inlet . . . . .	12
orbi_define_block_for_flow_injection . . . . .	13
orbi_export_data_to_excel . . . . .	14
orbi_filter_files . . . . .	15
orbi_filter_flagged_data . . . . .	16
orbi_filter_isotopocules . . . . .	17
orbi_filter_isox . . . . .	17
orbi_filter_satellite_peaks . . . . .	18
orbi_filter_scan_intensity . . . . .	18
orbi_filter_weak_isotopocules . . . . .	19
orbi_find_isox . . . . .	19
orbi_find_raw . . . . .	20
orbi_flag_outliers . . . . .	20
orbi_flag_satellite_peaks . . . . .	21
orbi_flag_weak_isotopocules . . . . .	22
orbi_get_blocks_info . . . . .	23
orbi_get_data . . . . .	24
orbi_get_example_files . . . . .	25
orbi_get_problems . . . . .	26
orbi_get_settings . . . . .	26
orbi_identify_isotopocules . . . . .	27
orbi_options . . . . .	28
orbi_plot_isotopocule_coverage . . . . .	29
orbi_plot_raw_data . . . . .	30
orbi_plot_satellite_peaks . . . . .	32
orbi_plot_shot_noise . . . . .	33
orbi_plot_spectra . . . . .	34
orbi_read_isox . . . . .	35
orbi_read_raw . . . . .	36
orbi_segment_blocks . . . . .	37
orbi_set_settings . . . . .	38
orbi_simplify_isox . . . . .	39
orbi_start_aggregator . . . . .	39
orbi_summarize_results . . . . .	41

---

`isoorbi-package`*isoorbi: Process Orbitrap Isotopocule Data*

---

## Description

Read and process isotopocule data from an Orbitrap Isotope Solutions mass spectrometer. Citation: Kantnerova et al. (Nature Protocols, 2024).

## Details

Resources:

- Website for the isoorbi package: <https://isoorbi.isoverse.org>
- Package options: [orbi\\_options](#)

## Author(s)

**Maintainer:** Caj Neubauer <[caj.neubauer@colorado.edu](mailto:caj.neubauer@colorado.edu)> ([ORCID](#)) [copyright holder]

Authors:

- Sebastian Kopf <[sebastian.kopf@colorado.edu](mailto:sebastian.kopf@colorado.edu)> ([ORCID](#))
- Kristýna Kantnerová <[kristyna.kantnerova@colorado.edu](mailto:kristyna.kantnerova@colorado.edu)> ([ORCID](#))

## See Also

Useful links:

- <https://isoorbi.isoverse.org/>
- <https://github.com/isoverse/isoorbi>
- Report bugs at <https://github.com/isoverse/isoorbi/issues>

---

`orbi_add_blocks_to_plot`*Plot blocks background*

---

## Description

This function can be used to add colored background to a plot of dual-inlet data where different colors signify different data types (data, startup time, changeover time, unused). Note that this function only works with continuous and pseudo-log y axis, not with log y axes.

**Usage**

```

orbi_add_blocks_to_plot(
  plot,
  x = c("guess", "scan.no", "time.min"),
  data_only = FALSE,
  fill = .data$data_type,
  fill_colors = c("#1B9E77", "#D95F02", "#7570B3", "#E7298A", "#66A61E", "#E6AB02",
    "#A6761D", "#666666"),
  fill_scale = scale_fill_manual("blocks", values = fill_colors),
  alpha = 0.5,
  show.legend = !data_only
)

```

**Arguments**

plot	object with a dataset that has defined blocks
x	which x-axis to use (time vs. scan number). If set to "guess" (the default), the function will try to figure it out from the plot.
data_only	if set to TRUE, only the blocks flagged as "data" ( <code>orbi_get_option("data_type_data")</code> ) are highlighted
fill	what to use for the fill aesthetic, default is the block <code>data_type</code>
fill_colors	which colors to use, by default a color-blind friendly color palettes ( <code>RColorBrewer</code> , <code>dark2</code> )
fill_scale	use this parameter to replace the entire fill scale rather than just the <code>fill_colors</code>
alpha	opacity settings for the background
show.legend	whether to include the background information in the legend

---

orbi\_adjust\_block      *Manually adjust block delimiters*

---

**Description**

This function can be used to manually adjust where certain block starts or ends after it's been defined with `orbi_define_block_for_flow_injection()` or `orbi_define_blocks_for_dual_inlet()` using either time or scan number. Note that adjusting blocks removes all block segmentation. Make sure to call `orbi_segment_blocks()` **after** adjusting block delimiters.

**Usage**

```

orbi_adjust_block(
  dataset,
  block,
  filename = NULL,
  shift_start_time.min = NULL,

```

```

    shift_end_time.min = NULL,
    shift_start_scan.no = NULL,
    shift_end_scan.no = NULL,
    set_start_time.min = NULL,
    set_end_time.min = NULL,
    set_start_scan.no = NULL,
    set_end_scan.no = NULL
  )

```

## Arguments

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
block	the block for which to adjust the start and/or end
filename	needs to be specified only if the dataset has more than one filename
shift_start_time.min	if provided, the start time of the block will be shifted by this many minutes (use negative numbers to shift back)
shift_end_time.min	if provided, the end time of the block will be shifted by this many minutes (use negative numbers to shift back)
shift_start_scan.no	if provided, the start of the block will be shifted by this many scans (use negative numbers to shift back)
shift_end_scan.no	if provided, the end of the block will be shifted by this many scans (use negative numbers to shift back)
set_start_time.min	if provided, sets the start time of the block as close as possible to this time
set_end_time.min	if provided, sets the end time of the block as close as possible to this time
set_start_scan.no	if provided, sets the start of the block to this scan number (scan must exist in the dataset)
set_end_scan.no	if provided, sets the end of the block to this scan number (scan must exist in the dataset)

## Value

A data frame (tibble) with block limits altered according to the provided start/end change parameters. Any data that is no longer part of the original block will be marked with the value of `orbi_get_option("data_type_unused")`. Any previously applied segmentation will be discarded (segment column set to NA) to avoid unintended side effects.

---

orbi\_aggregate\_raw      *Aggregate data from raw files*

---

## Description

This function allows dynamic aggregation and validation of data read by `orbi_read_raw()`. Like `orbi_read_raw()`, it is designed to be fail safe by safely catching errors and reporting back on them (see `orbi_get_problems()`). This function should work out of the box for most files without additional modification of the aggregator.

## Usage

```
orbi_aggregate_raw(  
  files_data,  
  aggregator = "standard",  
  show_progress = rlang::is_interactive(),  
  show_problems = TRUE  
)
```

## Arguments

<code>files_data</code>	the files read in by <code>orbi_read_raw()</code>
<code>aggregator</code>	typically the name of a registered aggregator (see all with <code>orbi_get_option("aggregators")</code> ), default is the "standard" aggregator included in the package ( <code>orbi_get_aggregator("standard")</code> ). Other options are "minimal" ( <code>orbi_get_aggregator("minimal")</code> ) and "extended" ( <code>orbi_get_aggregator("extended")</code> ). The aggregator parameter can also directly be an aggregator tibble (created/modified with <code>orbi_start_aggregator()</code> and/or <code>orbi_add_to_aggregator()</code> ) that defines which data should be aggregated and how.
<code>show_progress</code>	whether to show a progress bar, by default always enabled when running interactively e.g. inside RStudio (and disabled in a notebook), turn off with <code>show_progress = FALSE</code>
<code>show_problems</code>	whether to show problems encountered along the way (rather than just keeping track of them with <code>orbi_get_problems()</code> ). Set to <code>show_problems = FALSE</code> to turn off the live printout. Either way, all encountered problems can be retrieved with running <code>orbi_get_problems()</code> for the returned list

## Value

a list of merged dataframes collected from the `files_data` based on the aggregator definitions

---

orbi\_analyze\_shot\_noise  
*Shot noise calculation*

---

### Description

This function computes the shot noise calculation.

### Usage

```
orbi_analyze_shot_noise(dataset, include_flagged_data = FALSE)
```

### Arguments

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
include_flagged_data	whether to include flagged data in the shot noise calculation (FALSE by default)

### Details

Analyze shot noise

will calculate for all combinations of filename, compound, and isotopocule in the provided dataset

### Value

The processed data frame with new columns: `n_effective_ions`, `ratio`, `ratio_rel_se.permil`, `shot_noise.permil`

---

orbi\_calculate\_ions     *Calculate ions from intensities*

---

### Description

This functions calculates ions (`ions.incremental`) from intensities based on the equation

$$N_{ions} = S/N \cdot C_N/z \cdot \sqrt{R_N/R} \cdot \sqrt{N_{MS}}$$

where S is the reported signal (intensity) of the isotopocule, N is the noise associated with the signal (peakNoise), measured at the resolution setting R (resolution), the noise factor  $C_N$  (CN) is the number of charges corresponding to the Orbitrap noise band at some reference resolution  $R_N$  (RN),  $N_{MS}$  is the number of microscans, and z is the charge per ion (charge) of the isotopocule. See Makarov and Denisov (2009) and Eiler et al. (2017) for details about this equation. The

default values for CN and RN are from the Orbitrap Exploris Isotope Solutions Getting Started Guide (BRE0032999, Revision A, October 2022). Note that the exact values of these factors are only critical if the number of ions are interpreted outside of ratio calculations (in ratio calculations, these factors cancel).

### Usage

```
orbi_calculate_ions(dataset, CN = 3, RN = 240000)
```

### Arguments

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
CN	noise factor
RN	reference resolution of the noise factor

### Details

If using a dataset read from isox files you might have to add a charge column if it does not yet exist that indicates the charge of the isotopocule. If using data from raw files, `orbi_identify_isotopocules()` will automatically call this function with the default CN and RN so you don't need to call it explicitly unless you want to change these parameters.

### Value

same object as provided in dataset with new column `ions.incremental`

---

`orbi_calculate_ratios` *Calculate direct isotopocule ratios*

---

### Description

This function calculates isotopocule/base peak ratios for all isotopocules. It does not summarize or average the ratios in any way. For a summarizing version of this function, see `orbi_summarize_results()`.

### Usage

```
orbi_calculate_ratios(dataset)
```

### Arguments

dataset	A data frame output after running <code>orbi_define_basepeak()</code>
---------	-----------------------------------------------------------------------

### Value

Returns a mutated dataset with `ratio` column added.

---

```
orbi_calculate_summarized_ratio
```

*Calculate isotopocule ratio*

---

### Description

This function calculates the ratio of two isotopocules (the numerator and denominator). This function averages multiple measurements of each using the `ratio_method` and returns a single value. Normally this function is not called directly by the user, but via the function `orbi_summarize_results()`, which calculates isotopocule ratios and other results for an entire dataset.

### Usage

```
orbi_calculate_summarized_ratio(  
  numerator,  
  denominator,  
  ratio_method = c("direct", "mean", "sum", "median", "geometric_mean", "slope",  
                  "weighted_sum")  
)
```

### Arguments

<code>numerator</code>	Column(s) used as numerator; contains ion counts
<code>denominator</code>	Column used as denominator; contains ion counts
<code>ratio_method</code>	Method for computing the ratio. <b>Please note well:</b> the formula used to calculate ion ratios matters! Do not simply use arithmetic mean. The best option may depend on the type of data you are processing (e.g., MS1 versus M+1 fragmentation). <code>ratio_method</code> can be one of the following: <ul style="list-style-type: none"><li>• <code>mean</code>: arithmetic mean of ratios from individual scans.</li><li>• <code>sum</code>: sum of all ions of the numerator across all scans divided by the sum of all ions observed for the denominator across all scans.</li><li>• <code>geometric_mean</code>: geometric mean of ratios from individual scans.</li><li>• <code>slope</code>: The ratio is calculated using the slope obtained from a linear regression model that is weighted by the numerator <code>x</code>, using stats: <code>:lm(x ~ y + 0, weights = x)</code>.</li><li>• <code>weighted_sum</code>: A derivative of the sum option. The weighing function ensures that each scan contributes equal weight to the ratio calculation, i.e. scans with more ions in the Orbitrap do not contribute disproportionately to the total sum of <code>x</code> and <code>y</code> that is used to calculate <code>x/y</code>.</li></ul>

### Value

Single value ratio between the isotopocules defined as numerator and denominator calculated using the `ratio_method`.

**Examples**

```
df <-
  system.file("extdata", "testfile_flow.isox", package = "isoorbi") |>
  orbi_read_isox()

ions_180 <- dplyr::filter(df, isotopocule == "180")$ions.incremental
ions_M0 <- dplyr::filter(df, isotopocule == "M0")$ions.incremental

orbi_calculate_summarized_ratio(
  numerator = ions_180, denominator = ions_M0, ratio_method = "sum"
)

orbi_calculate_summarized_ratio(
  numerator = ions_180, denominator = ions_M0, ratio_method = "slope"
)
```

---

orbi\_check\_isoraw      *Check for the isoorbi raw file reader*

---

**Description**

By default, this will install the isoraw reader if it is missing or outdated, and will ask the user to agree to Thermo's license agreement for the **Thermo RawFileReader** before proceeding. This function runs automatically during a raw file read and does not usually need to be called directly by the user.

**Usage**

```
orbi_check_isoraw(
  install_if_missing = !on_cran(),
  reinstall_if_outdated = !on_cran(),
  reinstall_always = FALSE,
  min_version = "0.2.2",
  source = paste0("https://github.com/isoverse/isoorbi/releases/download/isoraw-v",
    min_version),
  accept_license = FALSE,
  ...
)
```

**Arguments**

```
install_if_missing
      install the reader if it's missing

reinstall_if_outdated
      install the reader if it's outdated (i.e. not at least min_version)

reinstall_always
      whether to (re-)install no matter what
```

min_version	the minimum version number required
source	the URL (or local path) where to find the raw file reader, by default this is the latests release of the executables on github
accept_license	explicitly accept Thermo's license agreement (if this is FALSE and the license has not previously been accepted, you will be asked about it)
...	passed on to download.file if (re-) installing the reader

---

orbi\_default\_theme     *Default isoorbi plotting theme*

---

### Description

Default isoorbi plotting theme

### Usage

```
orbi_default_theme(text_size = 16, facet_text_size = 20)
```

### Arguments

text_size	a font size for text
facet_text_size	a font size for facet text

### Value

ggplot theme object

---

orbi\_define\_basepeak     *Define the denominator for ratio calculation*

---

### Description

orbi\_define\_basepeak() sets one isotopocule in the data frame as the base peak (ratio denominator) and calculates the instantaneous isotope ratios against it.

### Usage

```
orbi_define_basepeak(dataset, basepeak_def)
```

### Arguments

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
basepeak_def	The isotopocule that gets defined as base peak, i.e. the denominator to calculate ratios

**Value**

same object as provided in dataset without the rows of the basepeak isotopocule and instead three new columns called basepeak, basepeak\_ions, and ratio holding the basepeak information and the isotope ratios vs. the base peak

**Examples**

```
fpath <- system.file("extdata", "testfile_flow.isoX", package = "isoorbi")
df <- orbi_read_isoX(file = fpath) |>
  orbi_simplify_isoX() |>
  orbi_define_basepeak(basepeak_def = "M0")
```

---

orbi\_define\_blocks\_for\_dual\_inlet

*Binning raw data into blocks for dual inlet analyses*

---

**Description**

This function sorts out (bins) data into individual blocks of reference, sample, changeover time, and startup time.

**Usage**

```
orbi_define_blocks_for_dual_inlet(
  dataset,
  ref_block_time.min,
  change_over_time.min,
  sample_block_time.min = ref_block_time.min,
  startup_time.min = 0,
  ref_block_name = orbi_get_option("di_ref_name"),
  sample_block_name = orbi_get_option("di_sample_name")
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <a href="#">orbi_identify_isotopocules()</a> as well as with a tibble from <a href="#">orbi_get_data(peaks = everything())</a> or when reading from an IsoX file)
ref_block_time.min	time where the signal is stable when reference is analyzed
change_over_time.min	time where the signal is unstable after switching from reference to sample or back
sample_block_time.min	time where the signal is stable when sample is analyzed

```

startup_time.min
    initial time to stabilize spray
ref_block_name  the name of the reference being measured
sample_block_name
    the name of the sample being measured

```

## Value

A data frame (tibble) with block annotations in the form of the additional columns described below:

- `data_group` is an integer that numbers each data group (whether that's startup, a sample block, a segment, etc.) in each file sequentially to uniquely identify groups of data that belong together - this column is NOT static (i.e. functions like `orbi_adjust_block()` and `orbi_segment_blocks()` will lead to renumbering) and should be used purely for grouping purposes in calculations and visualization
- `block` is an integer counting the data blocks in each file (0 is the startup block)
- `sample_name` is the name of the material being measured as defined by the `ref_block_name` and `sample_block_name` parameters
- `segment` is an integer defines segments within individual blocks - this will be NA until the optional `orbi_segment_blocks()` is called
- `data_type` is a text value describing the type of data in each `data_group` - for a list of the main categories, call `orbi_get_options("data_type")`

---

```
orbi_define_block_for_flow_injection
```

*Define data block for flow injection*

---

## Description

Define a data block by either start and end time or start and end scan number. If you want to make segments in the blocks (optional), note that this function - manually defining blocks - removes all block segmentation. Make sure to call `orbi_segment_blocks()` **only after** finishing block definitions.

## Usage

```

orbi_define_block_for_flow_injection(
  dataset,
  start_time.min = NULL,
  end_time.min = NULL,
  start_scan.no = NULL,
  end_scan.no = NULL,
  sample_name = NULL
)

```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
start_time.min	set the start time of the block
end_time.min	set the end time of the block
start_scan.no	set the start scan of the block
end_scan.no	set the end scan of the block
sample_name	if provided, will be used as the <code>sample_name</code> for the block

**Value**

A data frame (tibble) with block definition added. Any data that is not part of a block will be marked with the value of `orbi_get_option("data_type_unused")`. Any previously applied segmentation will be discarded (segment column set to NA) to avoid unintended side effects.

---

orbi\_export\_data\_to\_excel

*Export data to excel*

---

**Description**

This functions exports the dataset into an Excel file. If the dataset is aggregated data, use the `include` parameter to decide which part of the data to export.

**Usage**

```
orbi_export_data_to_excel(
  dataset,
  file,
  dbl_digits = 7,
  int_format = "0",
  dbl_format = sprintf(sprintf("%%.%sf", dbl_digits), 0),
  include = c("file_info", "summary", "scans", "peaks", "problems"),
  show_progress = rlang::is_interactive()
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
file	file path to export the file - recursively creates the directory if non-existent
dbl_digits	how many digits to show for dbls (all are exported)

int_format	the excel formatting style for integers
dbl_format	the excel formatting style for doubles (created automatically from the dbl_digits parameter)
include	which tibbles to include if dataset is aggregated data. By default includes all but spectra
show_progress	whether to show a progress bar, by default always enabled when running interactively e.g. inside RStudio (and disabled in a notebook), turn off with show_progress = FALSE

**Value**

returns dataset invisibly for use in pipes

---

orbi\_filter\_files      *Basic generic data files filter*

---

**Description**

This is a basic filter function for file names, compounds and time ranges. For filtering isotopocules, this function calls `orbi_filter_isotopocules()` internally (as of isoorbi version 1.5.0 `orbi_filter_isotopocules()` can also be used directly instead of via this function). Default value for all parameters is NULL, i.e. no filter is applied.

**Usage**

```
orbi_filter_files(
  dataset,
  filenames = NULL,
  compounds = NULL,
  isotopocules = NULL,
  time_min = NULL,
  time_max = NULL
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
filenames	Vector of file names to keep, keeps all if set to NULL (the default)
compounds	Vector of compounds to keep, keeps all if set to NULL (the default)
isotopocules	Vector of isotopocules to keep, keeps all if set to NULL (the default)
time_min	Minimum retention time in minutes ( <code>time.min</code> ), no minimum if set to NULL (the default)
time_max	Maximum retention time in minutes ( <code>time.min</code> ), no maximum if set to NULL (the default)

**Value**

Filtered tibble

**Examples**

```
fpath <- system.file("extdata", "testfile_flow.isox", package = "isoorbi")
df <-
  orbi_read_isox(file = fpath) |>
  orbi_simplify_isox() |>
  orbi_filter_files(
    filenames = c("s3744"),
    compounds = "HSO4-",
    isotopocules = c("M0", "34S", "180")
  )
```

---

orbi\_filter\_flagged\_data

*Filter out flagged data*

---

**Description**

This function filters out data that have been previously flagged using functions `orbi_flag_satellite_peaks()`, `orbi_flag_weak_isotopocules()`, and/or `orbi_flag_outliers()`. Note that this function is no longer necessary to call explicitly as `orbi_analyze_shot_noise()` and `orbi_summarize_results()` automatically exclude flagged data.

**Usage**

```
orbi_filter_flagged_data(dataset)
```

**Arguments**

`dataset` a tibble with previously flagged data from `orbi_flag_satellite_peaks()`, `orbi_flag_weak_isotopocules()`, and/or `orbi_flag_outliers()`

**Value**

a dataset with the flagged data filtered out

---

orbi\_filter\_isotopocules  
*Filter isotopocules*

---

### Description

This function helps filter out missing isotopocules, unidentified peaks, or select for specific isotopocule. It can be called any time after `orbi_identify_isotopocules()` or after reading from an isox file. By default (i.e. if run without setting any parameters), it removes unidentified peaks and missing isotopocules and keeps all others.

### Usage

```
orbi_filter_isotopocules(  
  dataset,  
  isotopocules = c(),  
  keep_missing = FALSE,  
  keep_unidentified = FALSE  
)
```

### Arguments

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
isotopocules	if provided, only these isotopocules will be kept
keep_missing	whether to keep missing isotopocules in the peaks list (i.e. those that should be there but are not), default is not to keep them
keep_unidentified	whether to keep unidentified isotopocules in the peaks list (i.e. peaks that have not been identified as a specific isotopocule), default is not to keep them

### Value

the dataset but filtered for these isotopocules

---

orbi\_filter\_isox      *Filter isox files*

---

### Description

**[Deprecated]** `orbi_filter_isox()` was renamed `orbi_filter_files()` to incorporate its wider scope for filtering both isox and raw files data.

**Usage**

```
orbi_filter_isox(...)
```

**Arguments**

... arguments passed on to [orbi\\_filter\\_files\(\)](#)

---

orbi\_filter\_satellite\_peaks

*Function replaced by [orbi\\_flag\\_satellite\\_peaks\(\)](#)*

---

**Description**

Function replaced by [orbi\\_flag\\_satellite\\_peaks\(\)](#)

**Usage**

```
orbi_filter_satellite_peaks(...)
```

**Arguments**

... parameters passed on to the new function [orbi\\_flag\\_satellite\\_peaks\(\)](#).

---

orbi\_filter\_scan\_intensity

*Function replaced by [orbi\\_flag\\_outliers\(\)](#)*

---

**Description**

Function replaced by [orbi\\_flag\\_outliers\(\)](#)

**Usage**

```
orbi_filter_scan_intensity(..., outlier_percent)
```

**Arguments**

... parameters passed on to the new function [orbi\\_flag\\_outliers\(\)](#).

outlier\_percent

outlier\_percent needs to be between 0 and 10, flags extreme scans based on TIC x injection time (i.e., ion intensity)

---

`orbi_filter_weak_isotopocules`*Function replaced by orbi\_flag\_weak\_isotopocules()*

---

**Description**

Function replaced by `orbi_flag_weak_isotopocules()`

**Usage**

```
orbi_filter_weak_isotopocules(...)
```

**Arguments**

... parameters passed on to the new function `orbi_flag_weak_isotopocules()`.

---

`orbi_find_isox`*Find isox files*

---

**Description**

Finds all `.isox` files in a folder.

**Usage**

```
orbi_find_isox(folder, recursive = TRUE)
```

**Arguments**

`folder` path to a folder with isox files  
`recursive` whether to find files recursively

**Examples**

```
# all .isox files provided with the isoorbi package  
orbi_find_isox(system.file("extdata", package = "isoorbi"))
```

---

orbi_find_raw	<i>Find raw files</i>
---------------	-----------------------

---

**Description**

Finds all .raw files in a folder.

**Usage**

```
orbi_find_raw(folder, pattern = NULL, include_cache = TRUE, recursive = TRUE)
```

**Arguments**

folder	path to a folder with raw files
pattern	provide a name pattern to find only specific raw files
include_cache	whether to include .raw.cache.zip folders in the absence of the corresponding .raw file so that copies of the cache are read even in the absence of the original raw files
recursive	whether to find files recursively

**Examples**

```
# all .raw files provided with the isoorbi package
orbi_find_raw(system.file("extdata", package = "isoorbi"))
```

---

orbi_flag_outliers	<i>Flag outlier scans</i>
--------------------	---------------------------

---

**Description**

This function flags outliers using one of the different methods provided by the parameters (to use multiple, please call this function several times sequentially). Note that this function evaluates outliers within each "uidx", "filename", and "injection" (for those of the columns that exist), and additionally within each "block" and "segment" if `by_block = TRUE`. In addition to any groupings already defined before calling this function using `dplyr`'s `group_by()`. It restores the original groupings in the returned dataset.

**Usage**

```
orbi_flag_outliers(
  dataset,
  agc_fold_cutoff = NA_real_,
  agc_window = c(),
  agc_absolute_cutoff = c(),
  by_block = TRUE
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
agc_fold_cutoff	flags scans with a fold cutoff based on the average number of ions in the Orbitrap analyzer. For example, <code>agc_fold_cutoff = 2</code> flags scans that have more than 2 times, or less than 1/2 times the average. TIC multiplied by injection time serves as an estimate for the number of ions in the Orbitrap.
agc_window	flags scans with a critically low or high number of ions in the Orbitrap analyzer. Provide a vector with 2 numbers <code>c(x, y)</code> flagging the lowest x percent and highest y percent. TIC multiplied by injection time serves as an estimate for the number of ions in the Orbitrap.
agc_absolute_cutoff	flags scans with a number of ions in the Orbitrap analyzer outside of an absolute range. Provide a vector with 2 numbers <code>c(x, y)</code> flagging data below x and above y of the TIC multiplied by injection time (which serves as an estimate for the number of ions in the Orbitrap).
by_block	if the dataset has block and segment definitions, should the outlier flag be evaluated within each block+segment or globally? default is within each block+segment, switch to globally by turning <code>by_block = FALSE</code>

**Value**

same object as provided in dataset with new columns `is_outlier` and `outlier_type` (if they don't already exist) that flags outliers identified by the method and provides the type of outlier (e.g. "2 fold agc cutoff"), respectively.

**Examples**

```
fpath <- system.file("extdata", "testfile_flow.isox", package = "isoorbi")
df <-
  orbi_read_isox(file = fpath) |>
  orbi_simplify_isox() |>
  orbi_flag_outliers(agc_window = c(1,99))
```

---

**orbi\_flag\_satellite\_peaks**

*Flag minor satellite peaks* Flags minor signals for an isotopocule that matches multiple peaks within its exact mass +/- tolerance interval in the same scan. These are often small satellite peaks generated by the Fourier transform. However, if there are satellite peaks of high intensity or very many satellite peaks, it can indicate that the *m/z* and tolerance setting used for identifying isotopocules need to be revisited. Visualize the flagged satellite peaks with `orbi_plot_satellite_peaks()`.

---

**Description**

Flag minor satellite peaks Flags minor signals for an isotopocule that matches multiple peaks within its exact mass +/- tolerance interval in the same scan. These are often small satellite peaks generated by the Fourier transform. However, if there are satellite peaks of high intensity or very many satellite peaks, it can indicate that the m/z and tolerance setting used for identifying isotopocules need to be revisited. Visualize the flagged satellite peaks with `orbi_plot_satellite_peaks()`.

**Usage**

```
orbi_flag_satellite_peaks(dataset)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
---------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Value**

same object as provided in dataset with new column `is_satellite_peak` that flags satellite peaks

**Examples**

```
fpath <- system.file("extdata", "testfile_flow.isox", package = "isoorbi")
df <-
  orbi_read_isox(file = fpath) |>
  orbi_simplify_isox() |>
  orbi_flag_satellite_peaks()
```

---

```
orbi_flag_weak_isotopocules
  Flag weak isotopocules
```

---

**Description**

This function flags isotopocules that are not detected in a minimum of `min_percent` of scans that then can be easily visualized with `orbi_plot_isotopocule_coverage()`. It evaluates weak isotopocules within each "uidx", "filename", "block", "segment" and "injection" (for those of the columns that exist), in addition to any groupings already defined before calling this function using dplyr's `group_by()`. It restores the original groupings in the returned data.

**Usage**

```
orbi_flag_weak_isotopocules(dataset, min_percent = 100)
```

## Arguments

dataset	A simplified IsoX data frame to be processed
min_percent	A number between 0 and 100 (inclusive). Isotopocule must be observed in at least this percentage of scans (please note: the percentage is defined relative to the most commonly observed isotopocule of each compound). The default is 100, the most stringent condition to ensure reliable isotopocule coverage and ratio calculations across data blocks. If you lower the default, be mindful of potential misinterpretations from using isotopocules that are very close to their detection limit within a datablock. For continuous flow operations it may be necessary to make data blocks smaller using <code>orbi_define_block_for_flow_injection()</code> and <code>orbi_adjust_block()</code> .

## Value

same object as provided in dataset with new column `is_weak_isotopocule` that flags weak isotopocules.

## Examples

```
fpath <- system.file("extdata", "testfile_flow.iso", package = "isoorbi")
df <- orbi_read_isoX(file = fpath) |>
  orbi_simplify_isoX() |>
  orbi_flag_weak_isotopocules(min_percent = 100)
```

---

orbi\_get\_blocks\_info *Summarize blocks info*

---

## Description

This function provides an overview table `blocks_info` which shows information on blocks in the dataset (block number, sample name, data type, scan number and start time where a block starts, and scan number and end time where a block ends).

## Usage

```
orbi_get_blocks_info(
  dataset,
  .by = c("uidx", "filename", "injection", "data_group", "block", "sample_name",
         "data_type", "segment")
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
.by	grouping columns for block info (akin to dplyr's .by parameter e.g. in <code>dplyr::summarize()</code> ). If not set by the user, all columns in the parameter's default values are used, if present in the dataset.

**Value**

a block summary or if no blocks defined yet, an empty tibble (with warning)

---

orbi_get_data	<i>Get data frame from aggregated data</i>
---------------	--------------------------------------------

---

**Description**

Retrieve a specific subset of the aggregated data into a single data frame by specifying which columns to take from each dataset (file\_info, scans, peaks, etc.) using `dplyr::select()` syntax. If data from more than one dataset is selected (e.g. some columns from scans AND some from peaks), the datasets are combined with an `dplyr::inner_join()` using the columns listed in by (only the ones actually in the datasets). Joins that would lead to duplicated data entries (i.e. many-to-many joins) are not allowed and will throw an error to avoid unexpected replications of individual datapoints. If you really want to do such a join, you'll have to do it manually.

**Usage**

```
orbi_get_data(
  aggregated_data,
  file_info = c("filename"),
  scans = NULL,
  peaks = NULL,
  spectra = NULL,
  problems = NULL,
  summary = NULL,
  by = c("uidx", "scan.no")
)
```

**Arguments**

aggregated_data	datasets aggregated from <code>orbi_aggregate_raw()</code>
file_info	columns to get from the aggregated file_info, all <code>dplyr::select()</code> syntax is supported
scans	columns to get from the aggregated scans, all <code>dplyr::select()</code> syntax is supported

peaks	columns to get from the aggregated peaks, all <code>dplyr::select()</code> syntax is supported
spectra	columns to get from the aggregated spectra, all <code>dplyr::select()</code> syntax is supported
problems	columns to get from the aggregated problems, all <code>dplyr::select()</code> syntax is supported
summary	columns to get from the summary calculated via <code>orbi_summarize_results()</code> , all <code>dplyr::select()</code> syntax is supported. Warning: it is not advisable to combine columns from summary with anything other than <code>file_info</code> as it will lead to duplicated datasets given that summary integrates across multiple scans.
by	which columns to look for when joining datasets together. Make sure to include the relevant by columns in the selections of the individual datasets so they are joined correctly. The default is usually sufficient

**Value**

a tibble

---

`orbi_get_example_files`

*Get available example file(s)*

---

**Description**

This function will provide the path(s) to example file(s). If a requested file is not yet available locally but is available on <https://github.com/isoverse/isodata>, it will download it from there into local storage. By default, it will download only cache files (.raw.cache.zip) instead of the original .raw files because the cache files are significantly smaller. To download the original raw files instead, use `download_raw_files = TRUE`.

**Usage**

```
orbi_get_example_files(
  filenames,
  download_raw_files = FALSE,
  download_always = FALSE
)
```

**Arguments**

`filenames` names of the example files

`download_raw_files` should the original raw files be downloaded? By default only cache files (raw.cache.zip) are downloaded as they are usually much smaller. However, they will not work for retrieving additional spectra. To download the original spectra, switch to `download_raw_files = TRUE`

`download_always` whether to download files anew even if they already exist locally

**Value**

file path(s) that can be passed directly to `orbi_read_raw()`

---

<code>orbi_get_problems</code>	<i>Retrieve parsing problems</i>
--------------------------------	----------------------------------

---

**Description**

This function retrieves parsing problems encountered during the reading and processing of files.

This function prints out parsing problems encountered during the reading and processing of files.

**Usage**

```
orbi_get_problems(obj, strip_ansi = TRUE)
```

```
orbi_show_problems(obj)
```

**Arguments**

`obj` data object that holds problems information

`strip_ansi` whether to remove ansi characters from the message, yes by default

**Value**

tibble data frame with a list of problems encountered during processing

---

<code>orbi_get_settings</code>	<i>Get all isoorbi package settings</i>
--------------------------------	-----------------------------------------

---

**Description****[Deprecated]**

`orbi_get_settings()` was renamed `orbi_get_options()` as part of `isoorbi` switching from 'settings' to 'options' to be consistent with base R naming conventions

**Usage**

```
orbi_get_settings(pattern = NULL)
```

**Arguments**

`pattern` passed on to `orbi_get_options()`

---

`orbi_identify_isotopocules`*Identify isotopocules*

---

## Description

Map the mass spectral peaks to specific isotopocules based on their mass.

## Usage

```
orbi_identify_isotopocules(  
  aggregated_data,  
  isotopocules,  
  default_tolerance = 1,  
  default_charge = 1  
)
```

## Arguments

<code>aggregated_data</code>	either data aggregated from <code>orbi_aggregate_raw()</code> or a straight-up tibble data frame of the peaks (e.g. retrieved via <code>orbi_get_data(peaks = everything())</code> ).
<code>isotopocules</code>	list of isotopocules to map, can be a data frame/tibble, a named vector such as <code>c("M0" = 61.9878, "15N" = 62.9850)</code> , or the name of a file to read from (.csv/.tsv/.xlsx are all supported). If provided as a tibble/file, the required columns are <code>isotopocule/isotopolog</code> and <code>mz/mass</code> (these alternative names for the columns, including uppercase versions, are recognized automatically). In addition, <code>tolerance/tolerance [mmu]/tolerance [mDa]</code> , <code>charge/z</code> , <code>#compound/compound</code> , and <code>fragment</code> are recognized, as well as any other (arbitrarily named) columns with additional information. Character columns in the <code>isotopocules</code> table (including <code>isotopocule</code> and <code>compound</code> ) are turned into factors with levels that preserve the order of isotopocules. That means that to change the order of isotopocules in downstream plotting functions, make sure to list them in the order you'd like them presented in. Note that if <code>tolerance/tolerance [mmu]/tolerance [mDa]</code> or <code>charge/z</code> are not provided, the values in the parameters <code>default_tolerance</code> and <code>default_charge</code> are used, respectively.
<code>default_tolerance</code>	tolerance (in mmu) to be used for isotopocule identification if a <code>tolerance/tolerance [mmu]/tolerance [mDa]</code> column is not included in <code>isotopocules</code>
<code>default_charge</code>	charge to be used for any unidentified peak, and if a <code>charge/z</code> column is not included in <code>isotopocules</code>

## Value

same object as provided in `aggregated_data` with added columns `compound` (if provided), `itc_uidx` (introduced unique isotopocule index), `isotopocule`, `mzExact`, `charge`, and `ions.incremental`

(via `orbi_calculate_ions()`), as well as any other additional information columns provided in isotopocules. Note that if the default CN and RN values of `orbi_calculate_ions()` are not the ones that should be used, simply run `orbi_calculate_ions()` explicitly afterwards. Also note that the information about columns that were NOT aggregated in previous steps is purposefully not preserved at this step.

---

orbi\_options

*Package options*


---

### Description

These options are best set via `orbi_options()` and queried via `orbi_get_option()`. However, the base functions `options()` and `getOption()` work as well but require an `isoorbi.` prefix (the package name and a dot) for the option name. Setting an option to a value of NULL means that the default is used. `orbi_get_options()` is available as an additional convenience function to retrieve a subset of options with a regular expression pattern.

### Usage

```
orbi_options(...)
```

```
orbi_get_options(pattern = NULL)
```

```
orbi_get_option(x)
```

### Arguments

...	set package options, syntax identical to <code>options()</code>
pattern	to retrieve multiple options (as a list) with a shared pattern
x	name of the specific option to retrieve

### Functions

- `orbi_options()`: set/get option values
- `orbi_get_options()`: get a subset of option values that fit a pattern
- `orbi_get_option()`: retrieve the current value of one option (option must be defined for the package)

### Options for the isoorbi package

- `di_ref_name`: the text label for dual inlet reference blocks
- `di_sample_name`: the text label for dual inlet sample blocks
- `data_type_data`: the text used to flag raw data as actually being data
- `data_type_startup`: the text used to flag raw data as being part of the startup
- `data_type_changeover`: the text used to flag raw data as being part of a changeover

- `data_type_unused`: the text used to flag raw data as being unused
- `aggregators`: data aggregators for pulling data out of raw files. The list of available aggregators is accessible via `orbi_get_option("aggregators")`. Individual aggregators are available via the shortcut helper function `orbi_get_aggregator("standard")`. Register new/overwrite existing aggregators via `orbi_register_aggregator()`.
- `debug`: turn on debug mode
- `auto_use_ansi`: whether to automatically enable correct rendering of stylized (ansi) output in HTML reports from notebooks that call `library(isoorbi)`. Can be turned off by calling `isoorbi::orbi_options(auto_use_ansi = FALSE)` **before** call `library(isoorbi)`.

### Examples

```
# All default options
orbi_get_options()

# Options that contain 'data' in the name
orbi_get_options("data")

# Specific option
orbi_get_option("data_type_unused")

# Change an option
orbi_options(data_type_unused = "flagged")
orbi_get_option("data_type_unused")

# Change back to default
orbi_options(data_type_unused = NULL)
orbi_get_option("data_type_unused")
```

---

```
orbi_plot_isotopocule_coverage
      Isotopocule coverage
```

---

### Description

The coverage of each isotopocule across scans/time is an important indicator for data completeness. These functions provide ways to summarize and visualize the isotopocule coverage in a dataset.

### Usage

```
orbi_plot_isotopocule_coverage(
  dataset,
  isotopocules = c(),
  x = c("scan.no", "time.min"),
  x_breaks = scales::breaks_pretty(5),
  add_data_blocks = TRUE
)
```

```
orbi_get_isotopocule_coverage(dataset)
```

### Arguments

dataset	a data frame or aggregated dataset with satellite peaks already identified (i.e. after <code>orbi_flag_satellite_peaks()</code> )
isotopocules	which isotopocules to visualize, if none provided will visualize all (this may take a long time or even crash your R session if there are too many isotopocules in the data set)
x	x-axis column for the plot, either "time.min" or "scan.no", default is "scan.no"
x_breaks	what breaks to use for the x axis, change to make more specifid tickmarks
add_data_blocks	add highlight for data blocks if there are any block definitions in the dataset (uses <code>orbi_add_blocks_to_plot()</code> ). To add blocks manually, set <code>add_data_blocks = FALSE</code> and manually call the <code>orbi_add_blocks_to_plot()</code> function afterwards.

### Value

a ggplot object  
summary data frame

### Functions

- `orbi_plot_isotopocule_coverage()`: visualizes isotope coverage. Weak isotopocules (if previously defined by `orbi_flag_weak_isotopocules()`) are highlighted in red.
- `orbi_get_isotopocule_coverage()`: calculates which stretches of the data have data for which isotopocules. This function is usually used indirectly by `orbi_plot_isotopocule_coverage()` but can be called directly to investigate isotopocule coverage.

---

orbi_plot_raw_data	<i>Visualize data</i>
--------------------	-----------------------

---

### Description

Call this function to visualize orbitrap data vs. time or scan number. The most common uses are `orbi_plot_raw_data(y = intensity)`, `orbi_plot_raw_data(y = ratio)`, and `orbi_plot_raw_data(y = tic * it.ms)`. If the selected y is peak-specific data (rather than scan-specific data like `tic * it.ms`), the `isotopocules` argument can be used to narrow down which isotopocules will be plotted. By default includes all isotopocules that have not been previously identified by `orbi_flag_weak_isotopocules()` (if already called on dataset).

**Usage**

```

orbi_plot_raw_data(
  dataset,
  isotopocules = c(),
  x = c("time.min", "scan.no"),
  x_breaks = scales::breaks_pretty(5),
  y,
  y_scale = c("raw", "linear", "pseudo-log", "log"),
  y_scale_sci_labels = TRUE,
  color = .data$isotopocule,
  colors = c("#1B9E77", "#D95F02", "#7570B3", "#E7298A", "#66A61E", "#E6AB02", "#A6761D",
    "#666666", "#BBBBBB"),
  color_scale = scale_color_manual(values = colors),
  add_data_blocks = TRUE,
  add_all_blocks = FALSE,
  show_outliers = TRUE
)

```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
isotopocules	which isotopocules to visualize, if none provided will visualize all (this may take a long time or even crash your R session if there are too many isotopocules in the data set)
x	x-axis column for the plot, either "time.min" or "scan.no", default is "scan.no"
x_breaks	what breaks to use for the x axis, change to make more specific tickmarks
y	expression for what to plot on the y-axis, e.g. <code>intensity</code> , <code>tic * it.ms</code> (pick one isotopocules as this is identical for different isotopocules), <code>ratio</code> . Depending on the variable, you may want to adjust the <code>y_scale</code> and potentially <code>y_scale_sci_labels</code> argument.
y_scale	what type of y scale to use: "log" scale, "pseudo-log" scale (smoothly transitions to linear scale around 0), "linear" scale, or "raw" (if you want to add a y scale to the plot manually instead)
y_scale_sci_labels	whether to render numbers with scientific exponential notation
color	expression for what to use for the color aesthetic, default is isotopocule
colors	which colors to use, by default a color-blind friendly color palettes (RColorBrewer, dark2)
color_scale	use this parameter to replace the entire color scale rather than just the colors
add_data_blocks	add highlight for data blocks if there are any block definitions in the dataset (uses <code>orbi_add_blocks_to_plot()</code> ). To add blocks manually, set <code>add_data_blocks = FALSE</code> and manually call the <code>orbi_add_blocks_to_plot()</code> function afterwards.

add\_all\_blocks add highlight for all blocks, not just data blocks (equivalent to the data\_only = FALSE argument in `orbi_add_blocks_to_plot()`)

show\_outliers whether to highlight data previously flagged as outliers by `orbi_flag_outliers()`

### Value

a ggplot object

---

orbi\_plot\_satellite\_peaks  
*Visualize satellite peaks*

---

### Description

Call this function any time after flagging the satellite peaks to see where they are. Use the `isotopocules` argument to focus on the specific isotopocules of interest.

### Usage

```
orbi_plot_satellite_peaks(
  dataset,
  isotopocules = c(),
  x = c("scan.no", "time.min"),
  y = c("ions.incremental", "intensity"),
  x_breaks = scales::breaks_pretty(5),
  y_scale = c("log", "pseudo-log", "linear", "raw"),
  y_scale_sci_labels = TRUE,
  colors = c("#1B9E77", "#D95F02", "#7570B3", "#E7298A", "#66A61E", "#E6AB02", "#A6761D",
    "#666666", "#BBBBBB"),
  color_scale = scale_color_manual(values = colors)
)
```

### Arguments

dataset a data frame or aggregated dataset with satellite peaks already identified (i.e. after `orbi_flag_satellite_peaks()`)

isotopocules which isotopocules to visualize, if none provided will visualize all (this may take a long time or even crash your R session if there are too many isotopocules in the data set)

x x-axis column for the plot, either "time.min" or "scan.no", default is "scan.no"

y y-axis column for the plot, typically either "ions.incremental" or "intensity", default is "ions.incremental" (falls back to "intensity" if "ions.incremental" has not been calculated yet for the provided dataset)

x\_breaks what breaks to use for the x axis, change to make more specific tickmarks

y_scale	what type of y scale to use: "log" scale, "pseudo-log" scale (smoothly transitions to linear scale around 0), "linear" scale, or "raw" (if you want to add a y scale to the plot manually instead)
y_scale_sci_labels	whether to render numbers with scientific exponential notation
colors	which colors to use, by default a color-blind friendly color palettes (RColorBrewer, dark2)
color_scale	use this parameter to replace the entire color scale rather than just the colors

**Value**

a ggplot object

---

orbi\_plot\_shot\_noise *Make a shot noise plot*

---

**Description**

This function creates a shot noise plot using a shotnoise data frame created by the [orbi\\_analyze\\_shot\\_noise\(\)](#) function.

**Usage**

```
orbi_plot_shot_noise(
  shotnoise,
  x = c("time.min", "n_effective_ions"),
  permil_target = NA_real_,
  color = "ratio_label",
  colors = c("#1B9E77", "#D95F02", "#7570B3", "#E7298A", "#66A61E", "#E6AB02", "#A6761D",
            "#666666")
)
```

**Arguments**

shotnoise	a shotnoise data frame
x	x-axis for the shot noise plot, either "time.min" or "n_effective_ions"
permil_target	highlight the target permil in the shotnoise plot
color	which column to use for the color aesthetic (must be a factor)
colors	which colors to use, by default a color-blind friendly color palettes (RColorBrewer, dark2)

**Details**

plot shot noise

**Value**

a ggplot object

---

orbi\_plot\_spectra      *Plot mass spectra*

---

### Description

This function visualizes mass spectra from aggregated raw file data. The spectra have to be previously read in with `include_spectra = c(1, 10, 100)` in `orbi_read_raw()`. By default, this function tries to visualize different isotopocule ranges (monoisotopic peak, M+1, M+2, M+3). To focus only on isotopocules of interest, run `orbi_identify_isotopocules()` and `orbi_filter_isotopocules()` first.

### Usage

```
orbi_plot_spectra(
  aggregated_data,
  mz_min = 0,
  mz_max = Inf,
  mz_base_peak = NULL,
  mz_focus_nominal_offsets = 0:4,
  max_scans = 6,
  max_files = 4,
  label_peaks = TRUE,
  show_filenames = TRUE,
  show_ref_and_lock_peaks = TRUE,
  show_focus_backgrounds = TRUE,
  background_colors = c("#1B9E77", "#D95F02", "#7570B3", "#E7298A", "#66A61E", "#E6AB02",
    "#A6761D", "#666666", "#BBBBBB")
)
```

### Arguments

<code>aggregated_data</code>	data aggregated by <code>orbi_aggregate_raw()</code> and, optionally, with isotopocules already identified by <code>orbi_identify_isotopocules()</code> , and (also optionally), already filtered with <code>orbi_filter_isotopocules()</code>
<code>mz_min</code>	which mz to start the main plot window at. By default include all.
<code>mz_max</code>	which mz to end the main plot window at. By default include all.
<code>mz_base_peak</code>	where is the base peak at (approximately)?. If not specified (the default), takes the largest peak in the <code>mz_min</code> to <code>mz_max</code> window.
<code>mz_focus_nominal_offsets</code>	which panels to visualize? 0 = whole spectrum, 1 = spectrum around monoisotopic peak + 1 mu (M+1), 2 = M+2, etc. By default includes the whole spectrum and up to +1, +2, +3, and +4 peaks (if they exist). To visualize only the whole spectrum, use <code>mz_focus_nominal_offsets = 0</code> . Likewise, to visualize only the area around the monoisotopic peak +1, provide <code>mz_focus_monimal_offsets = 1</code> (or = <code>c(1, 2)</code> for both +1 and +2 windows).

max_scans	spectra from how many scans to show at most. By default up to 6 (the number of available linetypes). To show only the spectrum from a single scan, set <code>max_scans = 1</code> . If more than 6 scan spectra are allowed (and there are more than 6 loaded in the <code>aggregated_data</code> ), turns of the linetype aesthetic.
max_files	spectra from how many files to show at most. Each file is shown as an additional line of panels.
label_peaks	whether to label the peaks in the M+1/2/3 panels. If isotopocules are already identified from <code>orbi_identify_isotopocules()</code> , uses the isotopocule names, otherwise the m/z values. Peaks that are missing (identified by <code>orbi_identify_isotopocules()</code> ) in all spectra are highlighted in red. To avoid labeling unidentified/missing peaks, run <code>orbi_filter_isotopocules()</code> first.
show_filenames	whether to show the filename in the first panel of each row (usually the full spectrum panel)
show_ref_and_lock_peaks	whether to show reference and lock mass peaks in the spectrum
show_focus_backgrounds	whether to highlight the M+x panels with specific background colors that match them with the mass bands highlighted in the first panel
background_colors	the colors to use for the background highlighting

---

orbi_read_isox	<i>Read IsoX file</i>
----------------	-----------------------

---

## Description

Read an IsoX dataput file (`.isox`) into a tibble data frame.

## Usage

```
orbi_read_isox(file)
```

## Arguments

`file` Path to the `.isox` file(s), single value or vector of paths

## Details

Additional information on the columns:

- `filename`: name of the original Thermo `.raw` file processed by IsoX
- `scan.no`: scan number
- `time.min`: acquisition or retention time in minutes
- `compound`: name of the compound (e.g., `NO3-`)
- `isotopocule`: name of the isotopocule (e.g., `15N`); called `isotopolog` in `.isox`

- `ions.incremental`: estimated number of ions, in increments since it is a calculated number
- `tic`: total ion current (TIC) of the scan
- `it.ms`: scan injection time (IT) in millisecond (ms)

### Value

A tibble containing at minimum the columns `filename`, `scan.no`, `time.min`, `compound`, `isotopocule`, `ions.incremental`, `tic`, `it.ms`

### Examples

```
fpath <- system.file("extdata", "testfile_dual_inlet.isox", package = "isoorbi")
df <- orbi_read_isox(file = fpath)
```

---

orbi_read_raw	<i>Read RAW files</i>
---------------	-----------------------

---

### Description

Read raw data files (`.raw`) from Orbitrap IRMS runs directly. This function extracts all available information and thus can be relatively slow (~1s / Mb on a typical personal computer) but with the caching this is only true the first time. The results can be used directly or, more typically, are aggregated with `orbi_aggregate_raw()` to safely extract the relevant information for downstream processing. This function is designed to be fail save by safely catching errors and reporting back on them (see `orbi_get_problems()`).

### Usage

```
orbi_read_raw(
  file_paths,
  show_progress = rlang::is_interactive(),
  show_problems = TRUE,
  include_spectra = FALSE,
  read_cache = TRUE,
  cache = TRUE,
  cache_spectra = cache,
  keep_cached_spectra = cache
)
```

### Arguments

<code>file_paths</code>	paths to the <code>.raw</code> file(s), single value or vector of paths. Use <code>orbi_find_raw()</code> to get all raw files in a folder.
<code>show_progress</code>	whether to show a progress bar, by default always enabled when running interactively e.g. inside RStudio (and disabled in a notebook), turn off with <code>show_progress = FALSE</code>

show_problems	whether to show problems encountered along the way (rather than just keeping track of them with <code>orbi_get_problems()</code> ). Set to <code>show_problems = FALSE</code> to turn off the live printout. Either way, all encountered problems can be retrieved with running <code>orbi_get_problems()</code> for the returned list
include_spectra	whether to include the spectral data from specific scans (e.g. <code>include_spectra = c(5, 100, 200)</code> reads out the spectra from scans 5, 100, and 200 for each file if they exist) or from all scans ( <code>include_spectra = TRUE</code> ). Including many or all scan spectra makes the read process slower (especially if <code>cache_spectra = FALSE</code> ) and the returned data frame tibble significantly larger. The default is <code>FALSE</code> (i.e. scan spectra are not returned).
read_cache	whether to read the file from cached .parquet files (if they exist) or anew
cache	whether to automatically cache the read raw files (writes highly efficient .parquet files in a folder with the same name as the file .cache appended)
cache_spectra	whether to automatically cache requested scan spectra (this can take up significant disc space), by default the same as <code>cache</code>
keep_cached_spectra	whether to keep the spectra from a raw file that were previously cached whenever <code>include_spectra</code> changes and requires reading the file anew. Having this <code>TRUE</code> (the default) makes it faster to iterate on code that changes which spectra to read but leads to larger cache files.

### Value

a tibble data frame where each row holds the file path and nested tibbles of datasets extracted from the raw file (typically `file_info`, `scans`, `peaks`, and `spectra`). This is the safest way to extract the data without needing to make assumptions about compatibility across files. Extract your data of interest from the tibble columns or use `orbi_aggregate_raw()` to extract safely across files.

---

orbi\_segment\_blocks     *Segment data blocks*

---

### Description

This step is optional and is intended to make it easy to explore the data within a sample or ref data block. Note that any raw data not identified with `data_type` set to "data" (`orbi_get_option("data_type_data")`) will stay unsegmented. This includes raw data flagged as "startup", "changeover", and "unused".

### Usage

```
orbi_segment_blocks(
  dataset,
  into_segments = NULL,
  by_scans = NULL,
  by_time_interval = NULL
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
into_segments	segment each data block into this many segments. The result will have exactly this number of segments for each data block except for if there are more segments requested than observations in a group (in which case each observation will be one segment)
by_scans	segment each data block into segments spanning this number of scans. The result will be approximately the requested number of scans per segment, depending on what is the most sensible distribution of the data. For example, in a hypothetical data block with 31 scans, if <code>by_scans = 10</code> , this function will create 3 segments with 11, 10 and 10 scans each (most evenly distributed), instead of 4 segments with 10, 10, 10, 1 (less evenly distributed).
by_time_interval	segment each data block into segments spanning this time interval. The result will have the requested time interval for all segments except usually the last one which is almost always shorter than the requested interval.

---

orbi\_set\_settings      *Set package settings*

---

**Description****[Deprecated]**

`orbi_set_settings()` was renamed `orbi_options()` as part of `isoorbi` switching from 'settings' to 'options' to be consistent with base R naming conventions

**Usage**

```
orbi_set_settings(...)
```

**Arguments**

...      named arguments to set specific options, passed on to `orbi_options()`

---

orbi\_simplify\_isox      *Simplify IsoX data*

---

### Description

Keep only columns that are directly relevant for isotopocule ratio analysis. This function is optional and does not affect any downstream function calls.

### Usage

```
orbi_simplify_isox(dataset, add = c())
```

### Arguments

dataset	IsoX data that is to be simplified
add	additional columns to keep

### Value

A tibble containing only the 9 columns: filepath, filename, scan.no, time.min, compound, isotopocule, ions.incremental, tic, it.ms, plus any additional columns defined in the add argument

### Examples

```
fpath <- system.file("extdata", "testfile_flow.isox", package="isoorbi")  
df <- orbi_read_isox(file = fpath) |> orbi_simplify_isox()
```

---

orbi\_start\_aggregator      *Dynamic data aggregator*

---

### Description

These functions allow definition of custom data aggregators for processing data extracted from raw files. An aggregator is run on each imported file and pulls together the relevant data users are interested in while making sure data formats are correct so that the aggregated data can be merged across several imported files for fast downstream processing.

**Usage**

```

orbi_start_aggregator(name)

orbi_add_to_aggregator(
  aggregator,
  dataset,
  column,
  source = column,
  default = NA,
  cast = "as.character",
  regexp = FALSE,
  func = NULL,
  args = NULL
)

orbi_register_aggregator(aggregator, name = attr(aggregator, "name"))

orbi_get_aggregator(name)

```

**Arguments**

name	a descriptive name for the aggregator. This name is automatically used as the default name when registering the aggregator via <code>orbi_register_aggregator()</code> .
aggregator	the aggregator table generated by <code>orbi_start_aggregator()</code> or passed from a previous call to <code>orbi_add_to_aggregator()</code> for constructing the entire aggregator by piping
dataset	the name of the dataset to aggregate from ( <code>file_info</code> , <code>scans</code> , <code>peaks</code> , <code>spectra</code> )
column	the name of the column in which data should be stored
source	single character column name or vector of column names (if alternatives could be the source) where in the dataset to find data for the column. If a vector of multiple column names is provided (e.g. <code>source = c("a1", "a2")</code> ), the first column name that's found during processing of a dataset will be used and passed to the function defined in <code>func</code> (if any) and then the one defined in <code>cast</code> . To provide multiple parameters from the data to <code>func</code> , define a list instead of a vector <code>source = list("a", "b", "c")</code> or if multiple alternative columns can be the source for any of the arguments, define as <code>source = list(c("a1", "a2"), "b", c("c1", "c2", "c3"))</code>
default	the default value if no source columns can be found or another error is encountered during aggregation. Note that the default value will also be processed with the function in <code>cast</code> to make sure it has the correct data type.
cast	what to cast the values of the resulting column to, most commonly <code>"as.character"</code> , <code>"as.integer"</code> , <code>"as.numeric"</code> , or <code>"as.factor"</code> . This is required to ensure all aggregated values have the correct data type.
regexp	whether source column names should be interpreted as a regular expressions for the purpose of finding the relevant column(s). Note if <code>regexp = TRUE</code> , the search for the source column always becomes case-insensitive so this can also

	be used for a direct match of a source column whose upper/lower casing can be unreliable. If a column is matched by a regexp and also by a direct aggregator rule, the direct aggregator rule takes precedence.
func	name of a processing function to apply before casting the value with the cast function. This is optional and can be used to conduct more elaborate preprocessing of a data or combining data from multiple source columns in the correct way (e.g. pasting together from multiple columns).
args	an optional list of arguments to pass to the func in addition to the values coming from the source column(s)

**Value**

an orbi aggregator tibble

**Functions**

- `orbi_start_aggregator()`: starts the aggregator
- `orbi_add_to_aggregator()`: add additional column to aggregate data for. Overwrites an existing aggregator entry for the same dataset and column if it already exists.
- `orbi_register_aggregator()`: register an aggregator in the isoorbi options so it can be retrieved with `orbi_get_aggregator()`
- `orbi_get_aggregator()`: retrieve a registered aggregator (get all aggregators with `orbi_get_option("aggregators")`)

---

orbi\_summarize\_results

*Generate the results table*

---

**Description**

Contains the logic to generate the results table. It passes the `ratio_method` parameter to the `orbi_calculate_summarized_ratio()` function for ratio calculations.

**Usage**

```
orbi_summarize_results(
  dataset,
  ratio_method = c("mean", "sum", "median", "geometric_mean", "slope", "weighted_sum"),
  .by = c("block", "sample_name", "segment", "data_group", "data_type", "injection"),
  include_flagged_data = FALSE,
  include_unused_data = FALSE
)
```

**Arguments**

dataset	An aggregated dataset or a data frame of peaks (i.e. works directly after <code>orbi_identify_isotopocules()</code> as well as with a tibble from <code>orbi_get_data(peaks = everything())</code> or when reading from an IsoX file)
ratio_method	Method for computing the ratio. <b>Please note well:</b> the formula used to calculate ion ratios matters! Do not simply use arithmetic mean. The best option may depend on the type of data you are processing (e.g., MS1 versus M+1 fragmentation). <code>ratio_method</code> can be one of the following: <ul style="list-style-type: none"> <li>• mean: arithmetic mean of ratios from individual scans.</li> <li>• sum: sum of all ions of the numerator across all scans divided by the sum of all ions observed for the denominator across all scans.</li> <li>• geometric_mean: geometric mean of ratios from individual scans.</li> <li>• slope: The ratio is calculated using the slope obtained from a linear regression model that is weighted by the numerator <code>x</code>, using <code>stats::lm(x ~ y + 0, weights = x)</code>.</li> <li>• weighted_sum: A derivative of the sum option. The weighing function ensures that each scan contributes equal weight to the ratio calculation, i.e. scans with more ions in the Orbitrap do not contribute disproportionately to the total sum of <code>x</code> and <code>y</code> that is used to calculate <code>x/y</code>.</li> </ul>
.by	additional grouping columns for the results summary (akin to <code>dplyr</code> 's <code>.by</code> parameter e.g. in <code>dplyr::summarize()</code> ). If not set by the user, all columns in the parameter's default values are used, if present in the dataset. Note that the order of these is also used to arrange the summary.
include_flagged_data	whether to include flagged data in the calculations (FALSE by default)
include_unused_data	whether to include unused data in the calculations (FALSE by default), in addition to peaks actually flagged as <code>orbi_get_option("data_type_data")</code>

**Value**

Returns a results summary table (as tibble if `dataset` is a tibble, as `dataset$tibble` if `dataset` is aggregated raw data) with the columns `filename`, `compound`, `isotopocule` and `basepeak` as well as the grouping columns from the `.by` parameter that are part of the input dataset. Additionally this function adds the following results columns: `start_scan.no`, `end_scan.no`, `start_time.min`, `mean_time.min`, `end_time.min`, `ratio`, `ratio_sem`, `ratio_relative_sem_permil`, `shot_noise_permil`, `No.of.Scans`, `minutes_to_1e6_ions`

- `ratio`: The isotope ratio between the isotopocule and the basepeak, calculated using the `ratio_method`
- `ratio_sem`: Standard error of the mean for the ratio
- `number_of_scans`: Number of scans used for the final ratio calculation
- `minutes_to_1e6_ions`: Time in minutes it would take to observe 1 million ions of the isotopocule used as numerator of the ratio calculation.
- `shot_noise_permil`: Estimate of the shot noise (more correctly thermal noise) of the reported ratio in permil.
- `ratio_relative_sem_permil`: Relative standard error of the reported ratio in permil

**Examples**

```
fpath <- system.file("extdata", "testfile_flow.isox", package = "isoorbi")
df <- orbi_read_isox(file = fpath) |>
  orbi_simplify_isox() |>
  orbi_define_basepeak("M0") |>
  orbi_summarize_results(ratio_method = "sum")
```

# Index

dplyr::inner\_join(), 24  
dplyr::select(), 24, 25  
dplyr::summarize(), 24, 42

get\_pkg\_options (orbi\_options), 28  
getOption(), 28

isoorbi (isoorbi-package), 3  
isoorbi-package, 3

options(), 28  
orbi\_add\_blocks\_to\_plot, 3  
orbi\_add\_blocks\_to\_plot(), 30–32  
orbi\_add\_to\_aggregator  
    (orbi\_start\_aggregator), 39  
orbi\_add\_to\_aggregator(), 6, 40  
orbi\_adjust\_block, 4  
orbi\_adjust\_block(), 13, 23  
orbi\_aggregate\_raw, 6  
orbi\_aggregate\_raw(), 24, 27, 34, 36, 37  
orbi\_analyze\_shot\_noise, 7  
orbi\_analyze\_shot\_noise(), 33  
orbi\_calculate\_ions, 7  
orbi\_calculate\_ions(), 28  
orbi\_calculate\_ratios, 8  
orbi\_calculate\_summarized\_ratio, 9  
orbi\_calculate\_summarized\_ratio(), 41  
orbi\_check\_isoraw, 10  
orbi\_default\_theme, 11  
orbi\_define\_basepeak, 11  
orbi\_define\_block\_for\_flow\_injection, 13  
orbi\_define\_block\_for\_flow\_injection(), 4, 23  
orbi\_define\_blocks\_for\_dual\_inlet, 12  
orbi\_define\_blocks\_for\_dual\_inlet(), 4  
orbi\_export\_data\_to\_excel, 14  
orbi\_filter\_files, 15  
orbi\_filter\_files(), 17, 18  
orbi\_filter\_flagged\_data, 16  
orbi\_filter\_isotopocules, 17  
orbi\_filter\_isotopocules(), 15, 34, 35  
orbi\_filter\_isox, 17  
orbi\_filter\_isox(), 17  
orbi\_filter\_satellite\_peaks, 18  
orbi\_filter\_scan\_intensity, 18  
orbi\_filter\_weak\_isotopocules, 19  
orbi\_find\_isox, 19  
orbi\_find\_raw, 20  
orbi\_find\_raw(), 36  
orbi\_flag\_outliers, 20  
orbi\_flag\_outliers(), 18, 32  
orbi\_flag\_satellite\_peaks, 21  
orbi\_flag\_satellite\_peaks(), 18, 30, 32  
orbi\_flag\_weak\_isotopocules, 22  
orbi\_flag\_weak\_isotopocules(), 30  
orbi\_get\_aggregator  
    (orbi\_start\_aggregator), 39  
orbi\_get\_aggregator(), 41  
orbi\_get\_blocks\_info, 23  
orbi\_get\_data, 24  
orbi\_get\_data(peaks = everything()), 5, 7, 8, 11, 12, 14, 15, 17, 21, 22, 24, 31, 38, 42  
orbi\_get\_example\_files, 25  
orbi\_get\_isotopocule\_coverage  
    (orbi\_plot\_isotopocule\_coverage), 29  
orbi\_get\_option (orbi\_options), 28  
orbi\_get\_option(), 28  
orbi\_get\_options (orbi\_options), 28  
orbi\_get\_options(), 26, 28  
orbi\_get\_problems, 26  
orbi\_get\_problems(), 6, 36, 37  
orbi\_get\_settings, 26  
orbi\_identify\_isotopocules, 27  
orbi\_identify\_isotopocules(), 5, 7, 8, 11, 12, 14, 15, 17, 21, 22, 24, 31, 34, 35, 38, 42

orbi\_options, [3](#), [28](#)  
orbi\_options(), [28](#), [38](#)  
orbi\_plot\_isotopocule\_coverage, [29](#)  
orbi\_plot\_isotopocule\_coverage(), [22](#),  
[30](#)  
orbi\_plot\_raw\_data, [30](#)  
orbi\_plot\_satellite\_peaks, [32](#)  
orbi\_plot\_satellite\_peaks(), [21](#), [22](#)  
orbi\_plot\_shot\_noise, [33](#)  
orbi\_plot\_spectra, [34](#)  
orbi\_read\_isox, [35](#)  
orbi\_read\_raw, [36](#)  
orbi\_read\_raw(), [6](#), [26](#), [34](#)  
orbi\_register\_aggregator  
    (orbi\_start\_aggregator), [39](#)  
orbi\_register\_aggregator(), [40](#)  
orbi\_segment\_blocks, [37](#)  
orbi\_segment\_blocks(), [4](#), [13](#)  
orbi\_set\_settings, [38](#)  
orbi\_show\_problems (orbi\_get\_problems),  
[26](#)  
orbi\_simplify\_isox, [39](#)  
orbi\_start\_aggregator, [39](#)  
orbi\_start\_aggregator(), [6](#), [40](#)  
orbi\_summarize\_results, [41](#)  
orbi\_summarize\_results(), [9](#), [25](#)