

# Package ‘tidyDenovix’

May 8, 2026

**Title** Cleans Spectrophotometry Data Obtained from the Denovix DS-11 Instrument

**Version** 2.1.0

**Description**

Cleans spectrophotometry data obtained from the Denovix instrument. The package also provides an option to normalize the data in order to compare the quality of the samples obtained.

**License** MIT + file LICENSE

**Encoding** UTF-8

**RoxygenNote** 7.2.3

**Suggests** knitr, learnr, rmarkdown, testthat (>= 3.0.0)

**Config/testthat/edition** 3

**Imports** data.table, janitor, readxl, tibble, dplyr, tidyr, stats, utils

**Depends** R (>= 2.10)

**LazyData** true

**URL** <https://github.com/AlphaPrime7/tidyDenovix>

**BugReports** <https://github.com/AlphaPrime7/tidyDenovix/issues>

**NeedsCompilation** no

**Author** Tingwei Adeck [aut, cre] (ORCID:  
<<https://orcid.org/0009-0009-7450-8863>>)

**Maintainer** Tingwei Adeck <[awesome.tingwei@outlook.com](mailto:awesome.tingwei@outlook.com)>

**Repository** CRAN

**Date/Publication** 2024-01-18 10:00:02 UTC

## Contents

extract_col_names	2
extract_sample_names	3
extract_wavelength	3

file_ext . . . . .	4
lambda_check . . . . .	5
lambda_check_source . . . . .	6
make_wavelength . . . . .	7
min_max_norm . . . . .	7
qc_attributes . . . . .	8
read_denovix . . . . .	9
read_denovix_data . . . . .	10
REnum . . . . .	10
rna1 . . . . .	11
rna2 . . . . .	11
rna3 . . . . .	12
spec . . . . .	12
tidyDenovix . . . . .	12

<b>Index</b>	<b>14</b>
--------------	-----------

---

extract_col_names	<i>Title: Extract key colnames from the Denovix data frame</i>
-------------------	--

---

## Description

Title: Extract key colnames from the Denovix data frame

## Usage

```
extract_col_names(xdf)
```

## Arguments

xdf                    The data frame for colname(s) extraction.

## Value

A vector of key column names.

## Author(s)

Tingwei Adeck

## Examples

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = read_denovix_data(fpath, file_type = 'csv')
col_names = extract_col_names(rna_data)
```

---

extract\_sample\_names *Title: Extract sample names from the Denovix data frame*

---

**Description**

Title: Extract sample names from the Denovix data frame

**Usage**

```
extract_sample_names(dfile, file_type = NULL)
```

**Arguments**

dfile            The denovix raw file for sample name(s) extraction.  
file\_type        The type of file.

**Value**

A vector of sample names.

**Author(s)**

Tingwei Adeck

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)  
esn = extract_sample_names(fpath, file_type = 'csv')
```

---

extract\_wavelength *Title: Extract wavelength*

---

**Description**

Title: Extract wavelength

**Usage**

```
extract_wavelength(xdf)
```

**Arguments**

xdf            The original data frame derived from importing Denovix data.

**Value**

A numeric data frame for the wavelength attribute.

**Author(s)**

Tingwei Adeck

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = read_denovix_data(fpath, file_type = 'csv')
wl = extract_wavelength(rna_data)
```

---

file\_ext

*Title: File Extension Finder*

---

**Description**

Title: File Extension Finder

**Usage**

```
file_ext(epath)
```

**Arguments**

epath            File path.

**Value**

A string representing the file extension.

**Author(s)**

Unknown (Adapted by Tingwei Adeck)

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
ext = file_ext(fpath)
```

---

lambda_check	<i>Title: Wavelength quality control</i>
--------------	--

---

## Description

Title: Wavelength quality control

## Usage

```
lambda_check(  
  qdf,  
  sample_type = c("RNA", "DNA"),  
  check_level = c("strict", "lax")  
)
```

## Arguments

qdf	A data frame with quality attributes.
sample_type	The type of sample under investigation.
check_level	The level of strictness based on sample type.

## Value

A data frame that meets the quality check criteria.

## Note

Some key assumptions are made about quality for RNA or DNA. At the moment column names is the main issue found with using this approach.

## Author(s)

Tingwei Adeck

## Examples

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)  
rna_data = read_denovix_data(fpath, file_type = 'csv')  
qc_check = lambda_check(rna_data, sample_type='RNA', check_level='lax')
```

lambda\_check\_source    *Title: Wavelength quality control*

---

## Description

Title: Wavelength quality control

## Usage

```
lambda_check_source(  
  odf,  
  sample_type = c("RNA", "DNA"),  
  check_level = c("strict", "lax")  
)
```

## Arguments

odf	A data frame with quality attributes.
sample_type	The type of sample under investigation.
check_level	The level of strictness based on sample type.

## Value

A vector of sample names for the different QC criteria.

## Note

Some key assumptions are made about quality for RNA or DNA. At the moment column names is the main issue found with using this approach.

## Author(s)

Tingwei Adeck

## Examples

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)  
rna_data = read_denovix_data(fpath, file_type = 'csv')  
qc_check = lambda_check_source(rna_data, sample_type='RNA', check_level='lax')
```

---

make_wavelength	<i>Title: Make wavelength</i>
-----------------	-------------------------------

---

**Description**

Title: Make wavelength

**Usage**

```
make_wavelength()
```

**Value**

A numeric data frame for the wavelength attribute.

**Author(s)**

Tingwei Adeck

**Examples**

```
w1 = make_wavelength()
```

---

min_max_norm	<i>Title: Min-Max normalization of attributes that require normalization</i>
--------------	--

---

**Description**

Title: Min-Max normalization of attributes that require normalization

**Usage**

```
min_max_norm(x)
```

**Arguments**

x                    A single value from an attribute passed in the function for normalization.

**Value**

A normalized value (value between 1 and 0)

**Note**

lapply is needed to apply the function across several columns in a data set.

**Author(s)**

Tingwei Adeck (Adapted from Statology)

**References**

<https://www.statology.org/how-to-normalize-data-in-r/>

**Examples**

```
test_df <- as.data.frame(c(seq(40)))
colnames(test_df) <- "test"
test_df_norm <- lapply(test_df[1:ncol(test_df)], min_max_norm)
```

---

qc\_attributes

*Title: Quality Control data frame*

---

**Description**

Title: Quality Control data frame

**Usage**

```
qc_attributes(dfile, file_type = NULL, xdf)
```

**Arguments**

dfile	The Denovix file path.
file_type	The type of file.
xdf	The Denovix data frame.

**Value**

A quality control data frame.

**Author(s)**

Tingwei Adeck

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = read_denovix_data(fpath, file_type = 'csv')
qc_attributes = qc_attributes(fpath, file_type = 'csv', rna_data)
```

---

read_denovix	<i>Title: Read Denovix files</i>
--------------	----------------------------------

---

**Description**

A function read Denovix data files.

**Usage**

```
read_denovix(dfile)
```

**Arguments**

dfile            A Denovix file or path to the Denovix file.

**Value**

A data frame.

**Note**

Denovix files can be saved as csv, txt or even excel files. This function accounts for these file types.

**Author(s)**

Tingwei Adeck

**See Also**

[read\\_denovix\\_data\(\)](#)

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = read_denovix(fpath)
```

read\_denovix\_data      *Title: Read Denovix files*

---

**Description**

A function read Denovix data files.

**Usage**

```
read_denovix_data(dfile, file_type = c("csv", "txt", "excel"))
```

**Arguments**

dfile                    A Denovix file or path to the Denovix file.  
file\_type                The file type being imported.

**Value**

A data frame.

**Note**

Denovix files can be saved as csv, txt or even excel files. This function accounts for these file types.

**Author(s)**

Tingwei Adeck

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = read_denovix_data(fpath, file_type = 'csv')
```

---

REnum                    *Title: Enum function*

---

**Description**

Title: Enum function

**Usage**

```
REnum(elist)
```

**Arguments**

elist                    A list with variable binding.

**Value**

The list supplied by the user.

**Author(s)**

Tingwei Adeck

**Examples**

```
cols = list('x260' = 1, 'x280' = 2, 'x230' = 3)
qc = REnum(cols)
a = (qc$x260)
```

---

rna1	<i>rna1.</i>
------	--------------

---

**Description**

An excel file version of the Denovix Spectrophotometry data.

**Usage**

```
rna1
```

**Format**

An object of class `tbl_df` (inherits from `tbl`, `data.frame`) with 14 rows and 151 columns.

---

rna2	<i>rna2.</i>
------	--------------

---

**Description**

A csv file version of the Denovix Spectrophotometry data.

**Usage**

```
rna2
```

**Format**

An object of class `data.frame` with 15 rows and 151 columns.

---

rna3	<i>rna3.</i>
------	--------------

---

**Description**

A csv file version of the Denovix Spectrophotometry data.

**Usage**

```
rna3
```

**Format**

An object of class `data.frame` with 15 rows and 152 columns.

---

spec	<i>spec.</i>
------	--------------

---

**Description**

A txt file version of the Denovix Spectrophotometry data.

**Usage**

```
spec
```

**Format**

An object of class `data.frame` with 15 rows and 151 columns.

---

tidyDenovix	<i>Title: Clean data from the Denovix DS-11 instrument</i>
-------------	--

---

**Description**

Title: Clean data from the Denovix DS-11 instrument

**Usage**

```
tidyDenovix(  
  dfile,  
  file_type = NULL,  
  sample_type = c("RNA", "DNA"),  
  check_level = c("strict", "lax"),  
  qc_omit = NULL,  
  normalized = c("yes", "no"),  
  fun = NA  
)
```

**Arguments**

<code>dfile</code>	The raw file obtained from the machine.
<code>file_type</code>	The file type specification.
<code>sample_type</code>	The sample type specification used in quality control.
<code>check_level</code>	The level of quality control performed.
<code>qc_omit</code>	Takes 'yes' or 'no' and determines if the qc data would be provided.
<code>normalized</code>	Takes 'yes' or 'no'.
<code>fun</code>	A parameter used for boolean expressions.

**Value**

A cleaned data frame with attribute names in some instances.

**Note**

The strict level of QC yields a data frame with no attribute names. This will be worked on so that users get the sample names for their data.

**Author(s)**

Tingwei Adeck

**Examples**

```
fpath <- system.file("extdata", "rnaspec2018.csv", package = "tidyDenovix", mustWork = TRUE)
rna_data = tidyDenovix(fpath, file_type = 'csv', sample_type = 'RNA', check_level = 'lax')
```

# Index

## \* datasets

rna1, [11](#)

rna2, [11](#)

rna3, [12](#)

spec, [12](#)

extract\_col\_names, [2](#)

extract\_sample\_names, [3](#)

extract\_wavelength, [3](#)

file\_ext, [4](#)

lambda\_check, [5](#)

lambda\_check\_source, [6](#)

make\_wavelength, [7](#)

min\_max\_norm, [7](#)

qc\_attributes, [8](#)

read\_denovix, [9](#)

read\_denovix\_data, [10](#)

read\_denovix\_data(), [9](#)

REnum, [10](#)

rna1, [11](#)

rna2, [11](#)

rna3, [12](#)

spec, [12](#)

tidyDenovix, [12](#)